

Coming
soon

MultiwellSlides & SlideFlask

Innovative slide systems for cell culture and microscopy

The new MultiwellSlides and SlideFlask family belong to a new product line available from August 2008. These microscope slide-based chamber systems are convenient to perform cell based applications (i.e. fluorescence staining) and cell culture testing directly on the slide.

Cultivating cells directly on the slide eliminates need for labour intensive transfer of cells from larger format tissue culture vessels to a slide format. Having a small growth area per well significantly reduces the amount of costly reagents for staining. MultiwellSlides with a lumox™, polyolefin and glass surface have standard microscope slide dimensions.


The polystyrene of the chamber of the lumox™ and polyolefin variant can be easily removed, which enables embedding slides with mounting media for staining conservation.

Applications include:

- Immunocytochemical assays
- Histological staining
- Cell culture assays

The MultiwellSlides are available in different formats (1-well, 2-well, 4-well and 8-wellSlide, SlideFlask) and different slide surfaces. The product portfolio encompasses the following:

- MultiwellSlides^{lumox™} and SlideFlask^{lumox™}
- MultiwellSlides^{polyolefin} and SlideFlask^{polyolefin}
- MultiwellSlides^{glass} and SlideFlask^{glass}
- MultiwellSlides^{coverglass I} and SlideFlask^{coverglass I}
- MultiwellSlides^{coverglass II} and SlideFlask^{coverglass II}
- MultiwellSlides^{polystyrene}



Features

- TC-treated for adherent cell culture
- Sterile, non-cytotoxic, DNase/RNase-free
- Slide materials include lumox™ surface, polyolefin, glass, coverglass and polystyrene
- Available in different well formats as well as in a SlideFlask format
- Optimally suited for cell culture testing, immunofluorescence applications and microscopy

Revision: March 2008 – 075 056

www.gbo.com/bioscience

Germany (Main office)
Greiner Bio-One GmbH
Phone: (+49) 7022 948-0
E-Mail: info@de.gbo.com

Belgium
Greiner Bio-One N. V.
Phone: (+32) 2-4 61 09 10
E-Mail: info@be.gbo.com

France
Greiner Bio-One SAS
Phone: (+33) 169-86 25 50
E-Mail: infos@fr.gbo.com

Japan
Greiner Bio-One Co. Ltd.
Phone: (+81) 3-35 05-88 75
E-Mail: info@jp.gbo.com

Austria
Greiner Bio-One GmbH
Phone: (+43) 7583 6791-0
E-Mail: office@at.gbo.com

Netherlands
Greiner Bio-One B. V.
Phone: (+31) 172-42 09 00
E-Mail: info@nl.gbo.com

UK
Greiner Bio-One Ltd.
Phone: (+44) 1453-82 52 55
E-Mail: info@uk.gbo.com

USA
Greiner Bio-One North America Inc.
Phone: (+1) 704-261-78 00
E-Mail: info@us.gbo.com



greiner bio-one

Order Information				
Cat.-No.	Product Description	Qty/ Bag	Qty/ Case	
9615 0001	SlideFlask ^{lumox™}	6	96	
9615 0101	1-wellSlide ^{lumox™}	6	96	
9615 0201	2-wellSlide ^{lumox™}	6	96	
9615 0401	4-wellSlide ^{lumox™}	6	96	
9615 0801	8-wellSlide ^{lumox™}	6	96	
9616 0101	1-wellSlide ^{polystyrene}	10	100	
9616 0201	2-wellSlide ^{polystyrene}	10	100	
9616 0401	4-wellSlide ^{polystyrene}	10	100	
9616 0801	8-wellSlide ^{polystyrene}	10	100	
9617 0002	SlideFlask ^{glass}	6	96	
9617 0102	1-wellSlide ^{glass}	6	96	
9617 0202	2-wellSlide ^{glass}	6	96	
9617 0402	4-wellSlide ^{glass}	6	96	
9617 0802	8-wellSlide ^{glass}	6	96	
9618 0002	SlideFlask ^{coverglass I}	6	96	
9618 0102	1-wellSlide ^{coverglass I}	6	96	
9618 0202	2-wellSlide ^{coverglass I}	6	96	
9618 0402	4-wellSlide ^{coverglass I}	6	96	
9618 0802	8-wellSlide ^{coverglass I}	6	96	
9619 0002	SlideFlask ^{coverglass II}	6	96	
9619 0102	1-wellSlide ^{coverglass II}	6	96	
9619 0202	2-wellSlide ^{coverglass II}	6	96	
9619 0402	4-wellSlide ^{coverglass II}	6	96	
9619 0802	8-wellSlide ^{coverglass II}	6	96	
9620 0001	SlideFlask ^{polyolefin}	6	96	
9620 0101	1-wellSlide ^{polyolefin}	6	96	
9620 0201	2-wellSlide ^{polyolefin}	6	96	
9620 0401	4-wellSlide ^{polyolefin}	6	96	
9620 0801	8-wellSlide ^{polyolefin}	6	96	

Revision: March 2008 – 075 056

www.gbo.com/bioscience

Germany (Main office)
Greiner Bio-One GmbH
Phone: (+49) 7022 948-0
E-Mail: info@de.gbo.com

Belgium
Greiner Bio-One N. V.
Phone: (+32) 2-4 61 09 10
E-Mail: info@be.gbo.com

France
Greiner Bio-One SAS
Phone: (+33) 169-86 25 50
E-Mail: infos@fr.gbo.com

Japan
Greiner Bio-One Co. Ltd.
Phone: (+81) 3-35 05-88 75
E-Mail: info@jp.gbo.com

Austria
Greiner Bio-One GmbH
Phone: (+43) 7583 6791-0
E-Mail: office@at.gbo.com

Netherlands
Greiner Bio-One B. V.
Phone: (+31) 172-42 09 00
E-Mail: info@nl.gbo.com

UK
Greiner Bio-One Ltd.
Phone: (+44) 1453-82 52 55
E-Mail: info@uk.gbo.com

USA
Greiner Bio-One North America Inc.
Phone: (+1) 704-261-78 00
E-Mail: info@us.gbo.com



greiner bio-one